

The Trust has supplied one surge protector for the 'KT Molecular Biology kit' and one for the 'KT Pathology kit'.

For *ALL* equipment requiring power and supplied by the Trust, you should endeavour to position it so that it can be plugged into the surge protector.

The following equipment **must** be used with the surge protector:

- Bio-Rad MyCycler thermal cycler
- Power pack
- Dell E5500 laptop
- UV Transilluminator
- Sigma 1-14 microcentrifuge
- pH meter
- Balance
- Microscopes – stereo and compound
- Microplate reader
- Spectrophotometer

